

AED/Bleed Kit Inspection Form

| | |
|--------------------|------------------------------------|
| Inspected By: | Location of Inspection (Building): |
| Date of Inspection | AED Serial Number: |

Use the following maintenance checklist when you check your unit.

Check the following Pass or Fail

| Check the following | PASS | FAIL | Comments |
|---|--------------------------|--------------------------|-----------------|
| Is the green check showing that the unit is ready to use? | <input type="checkbox"/> | <input type="checkbox"/> | |
| Is the unit clean, undamaged, and free of excessive wear? | <input type="checkbox"/> | <input type="checkbox"/> | |
| Are there any cracks or loose parts in the housing? | <input type="checkbox"/> | <input type="checkbox"/> | |
| Verify that electrodes are within their expiration date. | <input type="checkbox"/> | <input type="checkbox"/> | |
| Verify electrodes are connected to the unit and sealed in their package. Replace if expired. | <input type="checkbox"/> | <input type="checkbox"/> | Exp. Date: |